

Bariatric Quad Canes

HCPCS: E0105



FIGURE A



FIGURE B

- Lightweight, extra-strong base
- Manufactured with aluminum tubing, steel base
- Comfortable vinyl, contoured hand grips
- Four point base, combined with offset handle provides additional stability and support
- Aluminum handle, reinforced steel base
- Handle height adjusts from 30" to 39"
- Weight Capacity 500 lb

PRODUCTS

Description

10316-4	Quad Cane,HD,Sm.Base,Slvr.,4/cs
10317-4	Quad Cane,HD,LgBase,Slvr,4/cs

SPECIFICATIONS

Description

Warranty	Limited Lifetime
Handle Height	30" - 39"
Weight Capacity	500 lb