

## Aluminum Crutches with Accessories

HCPCS: E0114,E0114



FIGURE A



FIGURE B




FIGURE C

- Lighter than wood, the standard aluminum crutch provides stability and durability
- Double extruded center tube provides additional strength to weight-bearing area
- Pushpin adjustment makes it easy-to-adjust crutch height in 1" increments
- Non-skid, jumbo-size vinyl contoured tips provide excellent traction
- Easy wing nut adjustments for hand grips
- Comfortable, durable accessories are latex free
- Weight Capacity 350 lb on aluminum crutch (excluding pediatric)
- Piggy back label supplied on product hangtag (excluding pediatric 10416-1)

PRODUCTS	Description
10400-8	Crutch, Adult, 8pr/cs
10401-1	Crutch,Alum.Youth,1pr/bx
10401-8	Crutch, Youth, 8pr/cs
10402-8	Crutch, Tall, 8pr/cs
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Weight Capacity	350 lb