

# Avance® Y-Connector

## Negative Pressure Wound Therapy System

### Tubing



- Dual lumen
- Translucent

### Connectors

- Easy-to-use
- Compatible with new Avance® Tubing and Avance® ViewPad™

### Controlled Fluid Management Technology™ (CFM)

- Filter controlled air vent
- Simultaneous CFM<sup>1</sup>



### Avance® Y-Connector

- Y-Connector attaches up to two Avance dressings to one pump for large or multiple wounds
- Tube design with three attachment points
- CFM uses a controlled air intake to enable continuous movement of fluid from the wound bed to the canister<sup>2</sup>
- Simultaneous CFM for two wounds per Y-Connector, for increased functionality and safety
- Easy-to-use connectors for convenient attachment and disconnection

### Ordering Information

Apria Product Order Number	Part Number	Description	Pcs/Case
M056573	772009	Avance® Y-Connector (Sterile)	10

1. Data on file with Mölnlycke Healthcare.

2. Data on file with Mölnlycke Healthcare.

Kerlix is a trademark of Covidien.

Apria Healthcare LLC complies with applicable Federal civil rights laws and does not discriminate on the basis of race, color, national origin, age, disability, or sex.

Apria Healthcare LLC cumple con las leyes federales de derechos civiles aplicables y no discrimina por motivos de raza, color, nacionalidad, edad, discapacidad o sexo.

Apria Healthcare LLC 遵守適用的聯邦民權法律規定，不因種族、膚色、民族血統、年齡、殘障或性別而歧視任何人。

The Mölnlycke Health Care, Avance, CFM Technology, ViewPad, Mepitel, Mepiseal and Safetac trademarks, names and logo types are registered globally to one or more of the Mölnlycke Health Care Group of Companies. Distributed by Apria Healthcare LLC. © 2018 Mölnlycke Health Care AB. All rights reserved.

Visit us at [apria.com](http://apria.com)



APRIA HEALTHCARE®