

Avance® ViewPad™

Negative Pressure Wound Therapy System

Body

- Inspection window
- Strong adhesive
- Thin and flexible wings

Tubing

- Dual lumen
- Translucent



Connector

- Easy-to-use

Slide clamp

- Round edges
- Liquid-tight when closed

Avance® ViewPad™

- Dual lumen tube design: one air lumen provides a controlled air release that regulates flow in the second fluid lumen, which transports exudate
- Flexible wings enable sufficient adhesion and conformability for small wounds and areas difficult to dress¹
- Transparent inspection window for easy positioning and application
- Strong adhesive on ViewPad body reduces the risk of air leaks²

Other Features:

- Translucent tubing for monitoring fluid and exudate
- Controlled Fluid Management Technology™ (CFM) and dual lumen tubing allows an air intake that enables a continuous transfer of fluid from the wound bed to the canister³
- Clamp minimizes the risk of splashing when disconnecting tubes containing fluid or exudate

Ordering Information

Apria Product Order Number	Part Number	Description	Pcs/Case
M056561	662018	Avance® ViewPad™	10
M056567	664152	Avance® Foam Dressing Kit incl Film with Safetac® technology, ViewPad™ Small	5
M056568	664252	Avance® Foam Dressing Kit incl Film with Safetac® technology, ViewPad™ Medium	5
M056569	664352	Avance® Foam Dressing Kit incl Film with Safetac® technology, ViewPad™ Large	5
M056562	662152	Avance® Foam Dressing Kit incl Transparent Film, ViewPad™ Small	5
M056564	662252	Avance® Foam Dressing Kit incl Transparent Film, ViewPad™ Medium	5
M056565	662352	Avance® Foam Dressing Kit incl Transparent Film, ViewPad™ Large	5
M056566	662752	Avance® Foam Dressing Kit incl Transparent Film, ViewPad™ X-Large	5



Avance® Foam Dressing Kit with Safetac® technology



Avance® Foam Dressing Kit with Transparent Film

1. Data on file with Mölnlycke Health Care

2. Data on file with Mölnlycke Health Care

3. Data on file with Mölnlycke Health Care

Apria Healthcare LLC complies with applicable Federal civil rights laws and does not discriminate on the basis of race, color, national origin, age, disability, or sex.

Apria Healthcare LLC cumple con las leyes federales de derechos civiles aplicables y no discrimina por motivos de raza, color, nacionalidad, edad, discapacidad o sexo.

Apria Healthcare LLC 遵守適用的聯邦民權法律規定，不因種族、膚色、民族血統、年齡、殘障或性別而歧視任何人。

The Mölnlycke Health Care, Avance, CFM Technology, ViewPad, Mepitel, Mepiseal and Safetac trademarks, names and logo types are registered globally to one or more of the Mölnlycke Health Care Group of Companies. Distributed by Apria Healthcare LLC. © 2018 Mölnlycke Health Care AB. All rights reserved.

Visit us at apria.com



APRIA HEALTHCARE®